

# Endoscope Consumables

# 04



# Endoscope Consumables

## **table of contents**

04-1 Distal Caps / Hoods

04-2 Valves

04-3 Inlet Seals

04-4 Balloons

04-5 O-Rings

04-6 Others

04-7 VERSA Consumables

# PENTAX Medical

## Accessories Programme

### for caring for your PENTAX Medical endoscopes

#### Original is better

Finely tuned system components are essential for safe and successful endoscopy. That is why PENTAX Medical original accessories aim to provide the highest possible quality in every detail. The following pages describe our carefully compiled accessory programme to ensure safe and effective cleaning of your PENTAX Medical endoscope.



#### How to choose the correct accessories:

Before you decide on new accessories, please ensure that your PENTAX Medical endoscope is compatible with the desired accessory.

- The illustration below is intended to help you to identify the relevant areas of your endoscope quickly and simply. Our PENTAX Medical specialists will be pleased to advise you if you have any additional questions.
- All PENTAX Medical original accessories in this catalogue have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

#### Orientation

In order to make your choice of PENTAX Medical accessories as simple and reliable as possible, we have identified fundamental product features with simple-to-understand symbols displayed next to the relevant technical data.



**Autoclavable**



**For single use only**



**OE-A63, OE-A65 Single Use, Sterile Disposable Elevator Cap (DEC™)**



<b>Model Name</b> OE-A63	<b>Description</b> Single use sterile disposable elevator cap
Compatible Endoscopes	ED34-i10T2

<b>Model Name</b> OE-A65	<b>Description</b> Single use sterile disposable elevator cap
Compatible Endoscopes	ED32-i10



**OE-A55 Distal Rubber Tip (C.A.P.™)**



<b>Model Name</b> OE-A55	<b>Description</b> Distal rubber tip
Compatible Endoscopes	ED34-i10T



**OE-A58 Distal Rubber Hood**



<b>Model Name</b> OE-A58	<b>Description</b> Distal rubber hood
Compatible Endoscopes	EG-2990Zi



**OE-A59 Distal Rubber Hood**



<b>Model Name</b> OE-A59	<b>Description</b> Distal rubber hood
Compatible Endoscopes	EC-3890LZi    EC-3890FZi    EC-3890MZi



**OF-B205 Single Use Suction Control Valve**



Model Name	Description		
OF-B205	Suction control valve for the J10 series		
Compatible Endoscopes	EB11-J10	EB19-J10	VNL15-J10
	EB15-J10	EB19-J10U	



**OF-B179 Suction Control Valve**



Model Name	Description			
OF-B179	Suction control valve (red)			
Compatible Endoscopes	EB-1170K	EB-1575K	EB-1972K	ECY-1570K
	EB-1570AK	EB-1970AK	EB-1975K	ECY-1575K
	EB-1570K	EB-1970K	EB-1990i	VNL-1570STK
	EB-1572K	EB-1970TK	EB-1970UK	



**OF-B163 Irrigation Cessation Valve**



Model Name	Description	
OF-B163	Irrigation cessation valve (blue)	
Compatible Endoscopes	ECY-1570K	ECY-1575K



**OF-B140 Suction Control Valve**



Model Name	Description			
OF-B140	Suction control valve (red)			
Compatible Endoscopes	FB-15BS	FCY-15RBS	FI-10BS	FI-16RBS
	FB-15RBS	FHY-10RBS	FI-10RBS	FNL-15RP3
	FB-18BS	FHY-15RBS	FI-13BS	FUR-9RBS
	FB-18RBS	FI-9BS	FI-13RBS	
	FCP-9P	FI-9RBS	FI-16BS	

### OF-B53 Suction Control Valve



Model Name	Description			
OF-B53	Suction control valve			
Compatible Endoscopes	FB-15BS	FCY-15RBS	FI-10BS	FI-16RBS
	FB-15RBS	FHY-10RBS	FI-10RBS	FNL-15RP3
	FB-18BS	FHY-15RBS	FI-13BS	FUR-9RBS
	FB-18RBS	FI-9BS	FI-13RBS	
	FCP-9P	FI-9RBS	FI-16BS	



### OF-B136 Suction Control Valve



Model Name	Description		
OF-B136	Suction control valve (red)		
Compatible Endoscopes	FB-8V	FB-15V	FB-19TV
	FB-10V	FB-18V	FCN-15X



### OF-B137 Irrigation Cessation Valve



Model Name	Description
OF-B137	Irrigation cessation valve (blue)
Compatible Endoscopes	FCN-15X



### OF-B147 Irrigation activation valve

Model Name	Description
OF-B147	Irrigation activation valve
Compatible Endoscopes	FCN-15X

## OF-B120 Suction Control Valve



Model Name	Description			
OF-B120	Suction valve (red)			
Compatible Endoscopes	EC-2990Li	EC-3890LK	EC38-i10L	EG-2990i
	EC-2990Fi	EC-3890FK	EC38-i10F	EG-2990K
	EC-2990Mi	EC-3890MK	EC38-i10M	EG-2990Zi
	EC-3490Li	EC-3890FK2	EC38-i10F2	EG-3490K
	EC-3490Fi	EC-3890MK2	EC38-i10M2	EG-3890TK
	EC-3490Mi	EC-3890LZi	EC38-i10NL	EG16-K10
	EC-3490LK	EC-3890FZi	EC38-i10NF	EG17-J10
	EC-3490FK	EC-3890MZi	EC38-i10cL	EG27-i10
	EC-3490MK	EC-3890TLK	EC38-i10cF	EG29-i10
	EC-3490TLi	EC-3890TFK	EC38-i10cM	EG29-i10c
	EC-3490TFi	EC34-i10L	EC38-i10cF2	EG34-i10
	EC-3490TMi	EC34-i10F	ED-3490TK	FC-38LV
	EC-380LKp	EC34-i10M	ED32-i10	FC-38FV
	EC-380FKp	EC34-i10NF	ED34-i10T	FD-34V2
	EC-380FK2p	EC34-i10NL	ED34-i10T2	FG-16V
	EC-380MK2p	EC34-i10TL	EE17-J10	FG-24V
	EC-3890Li	EC34-i10TF	EG-1690K	FG-29V
	EC-3890Fi	EC34-i10TM	EG-2490K	FS-34V
	EC-3890Mi	EC34-i10cL	EG-2790i	VSB-2990i
	EC-3890Fi2	EC34-i10cF	EG-2790K	
	EC-3890Mi2	EC34-i10cM	EG-290Kp	

## OF-B188 Air/Water Feeding Valve



Model Name	Description			
OF-B188	Air/Water valve (purple)			
Compatible Endoscopes	EC-2990Li	EC-3890FK	EC34-i10cF	EG-1690K
	EC-2990Fi	EC-3890MK	EC34-i10cM	EG-2490K
	EC-2990Mi	EC-3890FK2	EC38-i10L	EG-2790i
	EC-3490Li	EC-3890MK2	EC38-i10F	EG-2790K
	EC-3490Fi	EC-3890LZi	EC38-i10M	EG-2990i
	EC-3490Mi	EC-3890FZi	EC38-i10F2	EG-2990K
	EC-3490LK	EC-3890MZi	EC38-i10M2	EG-2990Zi
	EC-3490FK	EC-3890TLK	EC38-i10NL	EG-3490K
	EC-3490MK	EC-3890TFK	EC38-i10NF	EG-3890TK
	EC-3490TLi	EC34-i10L	EC38-i10cL	EG16-K10
	EC-3490TFi	EC34-i10F	EC38-i10cF	EG17-J10
	EC-3490TMi	EC34-i10M	EC38-i10cM	EG27-i10
	EC-3890Li	EC34-i10NF	EC38-i10cF2	EG29-i10
	EC-3890Fi	EC34-i10NL	ED-3490TK	EG29-i10c
	EC-3890Mi	EC34-i10TL	ED32-i10	EG34-i10
	EC-3890Fi2	EC34-i10TF	ED34-i10T	VSB-2990i
	EC-3890Mi2	EC34-i10TM	ED34-i10T2	
	EC-3890LK	EC34-i10cL	EE17-J10	



### OF-B194 Gas/Water Feeding Valve



**Model Name**  
OF-B194

**Description**  
Gas/water feeding valve (purple)

Compatible Endoscopes	EC-2990Li	EC-3890LK	EC34-i10cF	EE17-J10
	EC-2990Fi	EC-3890FK	EC34-i10cM	EG-1690K
	EC-2990Mi	EC-3890MK	EC38-i10L	EG-2490K
	EC-3490Li	EC-3890FK2	EC38-i10F	EG-2790i
	EC-3490Fi	EC-3890MK2	EC38-i10M	EG-2790K
	EC-3490Mi	EC-3890LZi	EC38-i10F2	EG-2990i
	EC-3490LK	EC-3890FZi	EC38-i10M2	EG-2990K
	EC-3490FK	EC-3890MZi	EC38-i10NL	EG-2990Zi
	EC-3490MK	EC-3890TLK	EC38-i10NF	EG-3490K
	EC-3490TLi	EC-3890TFK	EC38-i10cL	EG-3890TK
	EC-3490TFi	EC34-i10L	EC38-i10cF	EG16-K10
	EC-3490TMi	EC34-i10F	EC38-i10cM	EG17-J10
	EC-3890Li	EC34-i10M	EC38-i10cF2	EG27-i10
	EC-3890Fi	EC34-i10TL	ED-3490TK	EG29-i10
	EC-3890Mi	EC34-i10TF	ED32-i10	EG29-i10c
	EC-3890Fi2	EC34-i10TM	ED34-i10T	EG34-i10
	EC-3890Mi2	EC34-i10cL	ED34-i10T2	VSB-2990i

### OF-B121 Air/Water Feeding Valve



**Model Name**  
OF-B121

**Description**  
Air/water valve (blue)

Compatible Endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	EG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V

### OF-B183 Air/Water Feeding Valve



**Model Name**  
OF-B183

**Description**  
Air/water valve (green)

Compatible Endoscopes	ES-3870K
--------------------------	----------

**OF-B177 Suction Control Valve**



<b>Model Name</b>	<b>Description</b>	
OF-B177	Suction valve (yellow)	
Compatible Endoscopes	EE-1580K	ES-3870K

**OF-B105 Irrigation Cessation Valve**



<b>Model Name</b>	<b>Description</b>		
OF-B105	Irrigation control valve		
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	

**OF-B151 Irrigation Cessation Valve**



<b>Model Name</b>	<b>Description</b>		
OF-B151	Irrigation (deactivation) valve		
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	

**OF-B152 Water Feeding Valve**



<b>Model Name</b>	<b>Description</b>		
OF-B152	Irrigation (activation) valve (blue)		
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	



### OF-B171 Suction Control Valve



Model Name	Description		
OF-B171	Suction control valve (yellow)		
Compatible Endoscopes	EG-3270UK	EG-3670URK	EG-3870UTK



### OF-B172 Air/Water/Balloon Feeding Valve



Model Name	Description		
OF-B172	Air/water/balloon feeding valve (green)		
Compatible Endoscopes	EG-3270UK	EG-3670URK	EG-3870UTK



### OF-B206 Suction/Balloon Suction Control Valve



Model Name	Description		
OF-B206	Suction/balloon suction control valve (red)		
Compatible Endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT



### OF-B207 Air/Water/Balloon Feeding Valve



Model Name	Description		
OF-B207	Air/water/balloon feeding valve (purple)		
Compatible Endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT


**OF-B215 PENTAX Medical PROFILE Single Use Endoscope Inlet Seal**


<b>Model Name</b>	<b>Description</b>			
OF-B215	Single use rubber inlet seal			
Compatible Endoscopes	EC-2990Li	EC34-i10M	EG16-K10	EG-3270UK
	EC-2990Fi	EC34-i10TL	EG17-J10	EG-3670URK
	EC-2990Mi	EC34-i10TF	EG27-i10	EG-3870UTK
	EC-3490Li	EC34-i10TM	EG29-i10	EG34-J10U
	EC-3490Fi	EC34-i10cL	EG29-i10c	EG36-J10UR
	EC-3490Mi	EC34-i10cF	EG34-i10	EG38-J10UT
	EC-3490LK	EC34-i10cM	FC-38LV	ES-3870K
	EC-3490FK	EC38-i10L	FC-38FV	FB-8V
	EC-3490MK	EC38-i10F	FD-34V2	FB-10V
	EC-3490TLi	EC38-i10M	FG-16V	FB-15V
	EC-3490TFi	EC38-i10F2	FG-24V	FB-18V
	EC-3490TMi	EC38-i10M2	FG-29V	FB-19TV
	EC-380LKp	EC38-i10NL	FS-34V	FB-15BS
	EC-380FKp	EC38-i10NF	VSB-2990i	FB-15RBS
	EC-380FK2p	EC38-i10cL	EB-1170K	FB-18BS
	EC-380MK2p	EC38-i10cF	EB-1570AK	FB-18RBS
	EC-3890Li	EC38-i10cM	EB-1570K	FCN-15X
	EC-3890Fi	EC38-i10cF2	EB-1572K	FCP-9P
	EC-3890Mi	ED-3490TK	EB-1575K	FCY-15RBS
	EC-3890Fi2	ED32-i10	EB-1970AK	FHY-10RBS
	EC-3890Mi2	ED34-i10T	EB-1970K	FHY-15RBS
	EC-3890LK	ED34-i10T2	EB-1970TK	FI-9BS
	EC-3890FK	EE17-J10	EB-1972K	FI-9RBS
	EC-3890MK	EG-1690K	EB-1975K	FI-10BS
	EC-3890FK2	EG-2490K	EB-1990i	FI-10RBS
	EC-3890MK2	EG-2790i	EB11-J10	FI-13BS
	EC-3890LZi	EG-2790K	EB15-J10	FI-13RBS
	EC-3890FZi	EG-290Kp	EB19-J10	FI-16BS
	EC-3890MZi	EG-2990i	EB-1970UK	FI-16RBS
	EC-3890TLK	EG-2990K	EB19-J10U	FNL-15RP3
	EC-3890TFK	EG-2990Zi	ECY-1570K	FUR-9RBS
	EC34-i10L	EG-3490K	ECY-1575K	VNL-1570STK
	EC34-i10F	EG-3890TK	EE-1580K	VNL15-J10

## OF-B190 Rubber Inlet Seal



**Model Name**  
OF-B190

**Description**  
Rubber inlet seal

Model Name	Description
OF-B190	Rubber inlet seal
Compatible Endoscopes	EC-2990Li EC34-i10NF EG16-K10 EG-3670URK
	EC-2990Fi EC34-i10NL EG17-J10 EG-3870UTK
	EC-2990Mi EC34-i10TL EG27-i10 EG34-J10U
	EC-3490Li EC34-i10TF EG29-i10 EG36-J10UR
	EC-3490Fi EC34-i10TM EG29-i10c EG38-J10UT
	EC-3490Mi EC34-i10cL EG34-i10 ES-3870K
	EC-3490LK EC34-i10cF FC-38LV FB-8V
	EC-3490FK EC34-i10cM FC-38FV FB-10V
	EC-3490MK EC38-i10L FD-34V2 FB-15V
	EC-3490TLi EC38-i10F FG-16V FB-18V
	EC-3490TFi EC38-i10M FG-24V FB-19TV
	EC-3490TMi EC38-i10F2 FG-29V FB-15BS
	EC-380LKp EC38-i10M2 FS-34V FB-15RBS
	EC-380FKp EC38-i10NL VSB-2990i FB-18BS
	EC-380FK2p EC38-i10NF EB-1170K FB-18RBS
	EC-380MK2p EC38-i10cL EB-1570AK FCN-15X
	EC-3890Li EC38-i10cF EB-1570K FCP-9P
	EC-3890Fi EC38-i10cM EB-1572K FCY-15RBS
	EC-3890Mi EC38-i10cF2 EB-1575K FHY-10RBS
	EC-3890Fi2 ED-3490TK EB-1970AK FHY-15RBS
	EC-3890Mi2 ED32-i10 EB-1970K FI-9BS
	EC-3890LK ED34-i10T EB-1970TK FI-9RBS
	EC-3890FK ED34-i10T2 EB-1972K FI-10BS
	EC-3890MK EE17-J10 EB-1975K FI-10RBS
	EC-3890FK2 EG-1690K EB-1990i FI-13BS
	EC-3890MK2 EG-2490K EB11-J10 FI-13RBS
	EC-3890LZi EG-2790i EB15-J10 FI-16BS
	EC-3890FZi EG-2790K EB19-J10 FI-16RBS
	EC-3890MZi EG-290Kp EB-1970UK FNL-15RP3
	EC-3890TLK EG-2990i EB19-J10U FUR-9RBS
	EC-3890TFK EG-2990K EGY-1570K VNL-1570STK
	EC34-i10L EG-2990Zi EGY-1575K VNL15-J10
	EC34-i10F EG-3490K EE-1580K
	EC34-i10M EG-3890TK EG-3270UK



**OE-A51 Balloons for PENTAX Ultrasound Endoscopes (20 pcs)**



<b>Model Name</b> OE-A51	<b>Description</b> Balloons for endoscopic ultrasound application	
<b>Compatible Endoscopes</b>	EG-3670URK	EG36-J10UR



**OE-A52 Balloons for PENTAX Ultrasound Endoscopes (20 pcs)**



<b>Model Name</b> OE-A52	<b>Description</b> Balloons for endoscopic ultrasound application	
<b>Compatible Endoscopes</b>	EB-1970UK	



**OE-A56 Balloons for PENTAX Ultrasound Endoscopes (20 pcs)**



<b>Model Name</b> OE-A56	<b>Description</b> Balloons for endoscopic ultrasound application	
<b>Compatible Endoscopes</b>	EG-3270UK	



**OE-A60 Single Use Balloon**



<b>Model Name</b> OE-A60	<b>Description</b> Balloons for endoscopic ultrasound application	
<b>Compatible Endoscopes</b>	EG34-J10U	EG38-J10UT



### OE-A61 PENTAX Medical Single Use Balloon



<b>Model Name</b> OE-A61	<b>Description</b> Single use balloon
Compatible Endoscopes	EB19-J10U



### OF-A67 Balloon-O-Ring Kit



<b>Model Name</b> OF-A67	<b>Description</b> Balloon and o-ring kit for EUS Gastroscope
Compatible Endoscopes	EG-3870UTK



### OE-A53 Balloon Installation Device



<b>Model Name</b> OE-A53	<b>Description</b> Balloon installation device
Compatible Endoscopes	EB-1970UK    EB19-J10U

**OF-B219 O-Ring Set**



<b>Model Name</b> OF-B219	<b>Description</b> O-ring set and installation device
Compatible Valves	OF-B207

**OF-B214 O-Ring Set**



<b>Model Name</b> OF-B214	<b>Description</b> A set of 3 pcs. of O-rings for the cleaning adapter
Compatible Valves	OF-B212

**OF-B204 O-Ring Set**

<b>Model Name</b> OF-B204	<b>Description</b> O-ring set
Compatible Valves	BS-LL1

**OF-B195 O Ring Set**

<b>Model Name</b> OF-B195	<b>Description</b> O-ring set
Compatible Valves	OF-B194      OE-B10

**OF-B192 O-Ring Set For Air/Water Feeding Valves**



<b>Model Name</b> OF-B192	<b>Description</b> O-ring set for air/water feeding valves
Compatible Valves	OF-B188

**OF-B191 O-Ring Set for Air/Water/Balloon Feeding Valves**



<b>Model Name</b> OF-B191	<b>Description</b> O-ring set for air/water/balloon feeding valves
Compatible Valves	OF-B172



### OF-B187 O-Ring Set for Suction Control Valves



<b>Model Name</b> OF-B187	<b>Description</b> To be attached to the suction/balloon suction control valve
Compatible Valves	OF-B171      OF-B206

### OF-B184 O-Ring Set For Air/Water Feeding Valves



<b>Model Name</b> OF-B184	<b>Description</b> O-Ring-set for A/W valve
Compatible Valves	OF-B183

### OF-B182 O-Ring Set

<b>Model Name</b> OF-B182	<b>Description</b> O-ring set
Compatible Valves	OF-B113

### OF-B181 O-Ring Set for Suction Control Valves



<b>Model Name</b> OF-B181	<b>Description</b> O-Ring-Set for suction valve, 3pcs
Compatible Valves	OF-B179

### OF-B143 O Ring Set

<b>Model Name</b> OF-B143	<b>Description</b> O-ring-set for A/W valve, 6 pcs
Compatible Valves	OF-B121

### OF-B127 O-Ring Set for Suction Control Valves



<b>Model Name</b> OF-B127	<b>Description</b> O-ring-set for suction valve, 3pcs	
Compatible Valves	OF-B120	OF-B177

### OF-B100 O-Ring Set For Suction Control/Irrigation Cessation Valves



<b>Model Name</b> OF-B100	<b>Description</b> O-ring-set for suction valve, 2 X 3 pcs		
Compatible Valves	OF-B140	OF-B151	OF-B152



### OE-C12 Water Jet Check Valve Adapter



Model Name	Description			
OE-C12	Forward water jet check valve adapter			
Compatible Endoscopes	EC-2990Li	EC-3890Fi2	EC34-i10NF	EC38-i10cL
	EC-2990Fi	EC-3890Mi2	EC34-i10NL	EC38-i10cF
	EC-2990Mi	EC-3890LK	EC34-i10TL	EC38-i10cM
	EC-3490Li	EC-3890FK	EC34-i10TF	EC38-i10cF2
	EC-3490Fi	EC-3890MK	EC34-i10TM	EG-2990i
	EC-3490Mi	EC-3890FK2	EC34-i10cL	EG-2990K
	EC-3490LK	EC-3890MK2	EC34-i10cF	EG-2990Zi
	EC-3490FK	EC-3890LZi	EC34-i10cM	EG-3490K
	EC-3490MK	EC-3890FZi	EC38-i10L	EG-3890TK
	EC-3490TLi	EC-3890MZi	EC38-i10F	EG29-i10
	EC-3490TFi	EC-3890TLK	EC38-i10M	EG29-i10c
	EC-3490TMi	EC-3890TFK	EC38-i10F2	EG34-i10
	EC-3890Li	EC34-i10L	EC38-i10M2	VSB-2990i
	EC-3890Fi	EC34-i10F	EC38-i10NL	ES-3870K
	EC-3890Mi	EC34-i10M	EC38-i10NF	



### OE-C14 Check Valve

Model Name	Description	
OE-C14	Rubber check valve (black)	
Compatible Valves	OE-C12	OF-B121



### OE-C15 Check Valve Set (10 pcs)



Model Name	Description	
OE-C15	Rubber check valve set (10 pcs) (black)	
Compatible Valves	OE-C12	OF-B121

## OF-B118 Water Jet Connector Cap



**Model Name**  
OF-B118

**Description**  
Water jet connector cap

Compatible  
Endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10TL	EC38-i10cM
EC-2990Fi	EC-3890Mi2	EC34-i10TF	EC38-i10cF2
EC-2990Mi	EC-3890LK	EC34-i10TM	EG-2990i
EC-3490Li	EC-3890FK	EC34-i10cL	EG-2990K
EC-3490Fi	EC-3890MK	EC34-i10cF	EG-2990Zi
EC-3490Mi	EC-3890FK2	EC34-i10cM	EG-3490K
EC-3490LK	EC-3890MK2	EC38-i10L	EG-3890TK
EC-3490FK	EC-3890LZi	EC38-i10F	EG29-i10
EC-3490MK	EC-3890FZi	EC38-i10M	EG29-i10c
EC-3490TLi	EC-3890MZi	EC38-i10F2	EG34-i10
EC-3490TFi	EC-3890TLK	EC38-i10M2	VSB-2990i
EC-3490TMi	EC-3890TFK	EC38-i10NL	ES-3870K
EC-3890Li	EC34-i10L	EC38-i10NF	
EC-3890Fi	EC34-i10F	EC38-i10cL	
EC-3890Mi	EC34-i10M	EC38-i10cF	

## OF-B113 Irrigation Tube



**Model Name**  
OF-B113

**Description**  
Irrigation tube for WJ channel

Compatible  
Endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10NF	EC38-i10cL
EC-2990Fi	EC-3890Mi2	EC34-i10NL	EC38-i10cF
EC-2990Mi	EC-3890LK	EC34-i10TL	EC38-i10cM
EC-3490Li	EC-3890FK	EC34-i10TF	EC38-i10cF2
EC-3490Fi	EC-3890MK	EC34-i10TM	EG-2990i
EC-3490Mi	EC-3890FK2	EC34-i10cL	EG-2990K
EC-3490LK	EC-3890MK2	EC34-i10cF	EG-2990Zi
EC-3490FK	EC-3890LZi	EC34-i10cM	EG-3490K
EC-3490MK	EC-3890FZi	EC38-i10L	EG-3890TK
EC-3490TLi	EC-3890MZi	EC38-i10F	EG29-i10
EC-3490TFi	EC-3890TLK	EC38-i10M	EG29-i10c
EC-3490TMi	EC-3890TFK	EC38-i10F2	EG34-i10
EC-3890Li	EC34-i10L	EC38-i10M2	VSB-2990i
EC-3890Fi	EC34-i10F	EC38-i10NL	ES-3870K
EC-3890Mi	EC34-i10M	EC38-i10NF	

## OF-G11 Gas Adapter

Model Name	Description			
OF-G11	CO <sub>2</sub> gas adapter			
Compatible Endoscopes	EC-2990Li	EC-3890LK	EC38-i10L	EG-2990i
	EC-2990Fi	EC-3890FK	EC38-i10F	EG-2990K
	EC-2990Mi	EC-3890MK	EC38-i10M	EG-2990Zi
	EC-3490Li	EC-3890FK2	EC38-i10F2	EG-3490K
	EC-3490Fi	EC-3890MK2	EC38-i10M2	EG-3890TK
	EC-3490Mi	EC-3890LZi	EC38-i10NL	EG16-K10
	EC-3490LK	EC-3890FZi	EC38-i10NF	EG17-J10
	EC-3490FK	EC-3890MZi	EC38-i10cL	EG27-i10
	EC-3490MK	EC-3890TLK	EC38-i10cF	EG29-i10
	EC-3490TLi	EC-3890TFK	EC38-i10cM	EG29-i10c
	EC-3490TFi	EC34-i10L	EC38-i10cF2	EG34-i10
	EC-3490TMi	EC34-i10F	ED-3490TK	VSB-2990i
	EC-380LKp	EC34-i10M	ED32-i10	EE-1580K
	EC-380FKp	EC34-i10NF	ED34-i10T	EG-3270UK
	EC-380FK2p	EC34-i10NL	ED34-i10T2	EG-3670URK
	EC-380MK2p	EC34-i10TL	EE17-J10	EG-3870UTK
	EC-3890Li	EC34-i10TF	EG-1690K	EG34-J10U
	EC-3890Fi	EC34-i10TM	EG-2490K	EG36-J10UR
	EC-3890Mi	EC34-i10cL	EG-2790i	EG38-J10UT
	EC-3890Fi2	EC34-i10cF	EG-2790K	ES-3870K
EC-3890Mi2	EC34-i10cM	EG-290Kp		



## OE-B5 Injection Tube



Model Name	Description			
OE-B5	Injection tube for ultrasound endoscopes			
Compatible Endoscopes	EB-1970UK	EB19-J10U	EG-3270UK	EG-3870UTK



## OE-B3 Elevator Wire Cleaning Inlet Cap



Model Name	Description			
OE-B3	Elevator wire cleaning inlet cap			
Compatible Endoscopes	EB-1970UK	EB19-J10U	EG-3270UK	EG-3870UTK

## OF-G1 PENTAX Standard Adapter Sleeve



Model Name	Description			
OF-G1	Adapter sleeve			
Compatible Endoscopes	FC-38LV	FB-8V	FNL-7RP3	VNL-1590STi
	FC-38FV	FB-10V	FNL-10RP3	VNL8-J10
	FD-34V2	FB-15V	FNL-15RP3	VNL11-J10
	FG-16V	FB-18V	FUR-9RBS	VNL15-J10
	FG-24V	FB-19TV	VNL-1070STK	BS-LC1
	FG-29V	FCN-15X	VNL-1190STK	
	FS-34V	FCP-9P	VNL-1570STK	

## OF-B55 T-adapter

Model Name	Description			
OF-B55	T-adapter for luer-lock			
Compatible Endoscopes	EC-2990Li	EC34-i10NF	EG16-K10	EG-3670URK
	EC-2990Fi	EC34-i10NL	EG17-J10	EG-3870UTK
	EC-2990Mi	EC34-i10TL	EG27-i10	EG34-J10U
	EC-3490Li	EC34-i10TF	EG29-i10	EG36-J10UR
	EC-3490Fi	EC34-i10TM	EG29-i10c	EG38-J10UT
	EC-3490Mi	EC34-i10cL	EG34-i10	ES-3870K
	EC-3490LK	EC34-i10cF	FC-38LV	FB-8V
	EC-3490FK	EC34-i10cM	FC-38FV	FB-10V
	EC-3490MK	EC38-i10L	FD-34V2	FB-15V
	EC-3490TLi	EC38-i10F	FG-16V	FB-18V
	EC-3490TFi	EC38-i10M	FG-24V	FB-19TV
	EC-3490TMi	EC38-i10F2	FG-29V	FB-15BS
	EC-380LKp	EC38-i10M2	FS-34V	FB-15RBS
	EC-380FKp	EC38-i10NL	VSB-2990i	FB-18BS
	EC-380FK2p	EC38-i10NF	EB-1170K	FB-18RBS
	EC-380MK2p	EC38-i10cL	EB-1570AK	FCN-15X
	EC-3890Li	EC38-i10cF	EB-1570K	FCP-9P
	EC-3890Fi	EC38-i10cM	EB-1572K	FCY-15RBS
	EC-3890Mi	EC38-i10cF2	EB-1575K	FHY-10RBS
	EC-3890Fi2	ED-3490TK	EB-1970AK	FHY-15RBS
	EC-3890Mi2	ED32-i10	EB-1970K	FI-9BS
	EC-3890LK	ED34-i10T	EB-1970TK	FI-9RBS
	EC-3890FK	ED34-i10T2	EB-1972K	FI-10BS
	EC-3890MK	EE17-J10	EB-1975K	FI-10RBS
	EC-3890FK2	EG-1690K	EB-1990i	FI-13BS
	EC-3890MK2	EG-2490K	EB11-J10	FI-13RBS
	EC-3890LZi	EG-2790i	EB15-J10	FI-16BS
	EC-3890FZi	EG-2790K	EB19-J10	FI-16RBS
	EC-3890MZi	EG-290Kp	EB-1970UK	FNL-15RP3
	EC-3890TLK	EG-2990i	EB19-J10U	FUR-9RBS
	EC-3890TFK	EG-2990K	ECY-1570K	VNL-1570STK
	EC34-i10L	EG-2990Zi	ECY-1575K	VNL15-J10
	EC34-i10F	EG-3490K	EE-1580K	
EC34-i10M	EG-3890TK	EG-3270UK		



### OF-B216 Single use Irrigation Tap

Model Name	Description			
OF-B216	Single use irrigation tap (25pcs)			
Compatible Endoscopes	ECY-1570K	ECY-1575K	EG16-K10	EG-1690K



### OF-B161 Instrument Channel Selector



Model Name	Description		
OF-B161	Instrument channel selector for dual channel scopes		
Compatible Endoscopes	EC-3890TLK	EC-3890TFK	EG-3890TK

### OE-C28 Ventilation Cap



Model Name	Description			
OE-C28	Ventilation cap for IMAGINA endoscopes			
Compatible Endoscopes	EC34-i10cL	EC34-i10cM	EC38-i10cF	EC38-i10cF2
	EC34-i10cF	EC38-i10cL	EC38-i10cM	EG29-i10c

### OE-C29 Venting Connector Adapter



Model Name	Description			
OE-C29	Venting connector adapter for IMAGINA endoscopes			
Compatible Endoscopes	EC34-i10cL	EC34-i10cM	EC38-i10cF	EC38-i10cF2
	EC34-i10cF	EC38-i10cL	EC38-i10cM	EG29-i10c

### OE-C21 Venting Cap-English



Model Name	Description			
OE-C21	Venting cap-English			
Compatible Endoscopes	EC34-i10NF	EC34-i10TM	ED34-i10T	EG34-i10
	EC34-i10NL	EC38-i10NL	ED34-i10T2	
	EC34-i10TL	EC38-i10NF	EE17-J10	
	EC34-i10TF	ED32-i10	EG17-J10	

## OF-C5 Ventilation Cap



Model Name	Description			
OF-C5	Ventilation cap			
Compatible Endoscopes	EC-2990Li	EC34-i10F	EB-1572K	FCN-15X
	EC-2990Fi	EC34-i10M	EB-1575K	FCP-9P
	EC-2990Mi	EC38-i10L	EB-1970AK	FCY-15RBS
	EC-3490Li	EC38-i10F	EB-1970K	FHY-10RBS
	EC-3490Fi	EC38-i10M	EB-1970TK	FHY-15RBS
	EC-3490Mi	EC38-i10F2	EB-1972K	FI-7BS
	EC-3490LK	EC38-i10M2	EB-1975K	FI-7RBS
	EC-3490FK	ED-3490TK	EB-1990i	FI-9BS
	EC-3490MK	EG-1690K	EB11-J10	FI-9RBS
	EC-3490TLi	EG-2490K	EB15-J10	FI-10BS
	EC-3490TFi	EG-2790i	EB19-J10	FI-10RBS
	EC-3490TMi	EG-2790K	EB-1970UK	FI-13BS
	EC-380LKp	EG-290Kp	EB19-J10U	FI-13RBS
	EC-380FKp	EG-2990i	ECY-1570K	FI-16BS
	EC-380FK2p	EG-2990K	ECY-1575K	FI-16RBS
	EC-380MK2p	EG-2990Zi	EE-1580K	FNL-7RP3
	EC-3890Li	EG-3490K	EG-3270UK	FNL-10RBS
	EC-3890Fi	EG-3890TK	EG-3670URK	FNL-10RP3
	EC-3890Mi	EG16-K10	EG-3870UTK	FNL-15RP3
	EC-3890Fi2	EG27-i10	EG34-J10U	FUR-9RBS
	EC-3890Mi2	EG29-i10	EG36-J10UR	VNL-1070STK
	EC-3890LK	FC-38LV	EG38-J10UT	VNL-1190STK
	EC-3890FK	FC-38FV	ES-3870K	VNL-1570STK
	EC-3890MK	FD-34V2	FB-8V	VNL-1590STi
	EC-3890FK2	FG-16V	FB-10V	VNL8-J10
	EC-3890MK2	FG-24V	FB-15V	VNL11-J10
	EC-3890LZi	FG-29V	FB-18V	VNL15-J10
	EC-3890FZi	FS-34V	FB-19TV	PVK-1070Z
	EC-3890MZi	VSB-2990i	FB-15BS	PVK-J10
	EC-3890TLK	EB-1170K	FB-15RBS	
	EC-3890TFK	EB-1570AK	FB-18BS	
	EC34-i10L	EB-1570K	FB-18RBS	

## BS-LC1 Light Guide Cable



Model Name	Description			
BS-LC1	Light guide cable for BS-endoscopes			
Compatible Endoscopes	FB-15BS	FHY-15RBS	FI-10RBS	FUR-9RBS
	FB-15RBS	FI-7BS	FI-13BS	
	FB-18BS	FI-7RBS	FI-13RBS	
	FB-18RBS	FI-9BS	FI-16BS	
	FCY-15RBS	FI-9RBS	FI-16RBS	
	FHY-10RBS	FI-10BS	FNL-10RBS	



### BS-LL1 LED Light Source Batt-Type



Model Name	Description			
BS-LL1	LED light source battery type			
Compatible Endoscopes	FB-15BS	FHY-15RBS	FI-10RBS	FUR-9RBS
	FB-15RBS	FI-7BS	FI-13BS	
	FB-18BS	FI-7RBS	FI-13RBS	
	FB-18RBS	FI-9BS	FI-16BS	
	FCY-15RBS	FI-9RBS	FI-16RBS	
	FHY-10RBS	FI-10BS	FNL-10RBS	



### OF-Z5 Bite Block



Model Name	Description
OF-Z5	Bite block, autoclavable

### OF-Z11 Silicone Oil

Model Name	Description
OF-Z11	Silicone oil

### AP-PV1 Adapter for Video Module for PVK-1070Z



Model Name	Description			
AP-PV1	Adapter lens for PENTAX Medical video module			
Compatible Endoscopes	FC-38LV	FG-29V	FB-19TV	FI-7BS
	FC-38FV	FB-8V	FB-15BS	FI-9BS
	FD-34V2	FB-10V	FB-18BS	FI-10BS
	FG-16V	FB-15V	FCN-15X	FI-13BS
	FG-24V	FB-18V	FCP-9P	

### AP-RV3 Scope Mount Adapter



Model Name	Description
AP-RV3	Endoscope mount adapter for rigid ocular endoscopes
Compatible Endoscopes	PVK-J10

### AP-OV1 Scope Mount Adapter

Model Name	Description
AP-OV1	Adapter for Olympus OES scope to video module

### AP-RV1 Scope Mount Adapter for PVK-1070Z

<b>Model Name</b> AP-RV1	<b>Description</b> Adapter lens for PENTAX Medical video module		
Compatible Endoscopes	FB-15RBS	FB-18RBS	



### OF-B107 Protective Rubber Seal

<b>Model Name</b> OF-B107	<b>Description</b> Protective rubber seal		
Compatible Endoscopes	FCP-9P		

### AL-SL Storz Light Source Adapter

<b>Model Name</b> AL-SL	<b>Description</b> Storz light source adapter			
Compatible Endoscopes	FC-38LV	FB-8V	FNL-7RP3	VNL-1590STi
	FC-38FV	FB-10V	FNL-10RP3	VNL8-J10
	FD-34V2	FB-15V	FNL-15RP3	VNL11-J10
	FG-16V	FB-18V	FUR-9RBS	VNL15-J10
	FG-24V	FB-19TV	VNL-1070STK	BS-LC1
	FG-29V	FCN-15X	VNL-1190STK	
	FS-34V	FCP-9P	VNL-1570STK	

### AT-ST2 Ocular Adapter to Rigid Type

<b>Model Name</b> AT-ST2	<b>Description</b> Ocular adapter to rigid type			
Compatible Endoscopes	FC-38LV	FS-34V	FB-15BS	FI-10BS
	FC-38FV	FB-8V	FB-18BS	FI-13BS
	FD-34V2	FB-10V	FCN-15X	FI-16BS
	FG-16V	FB-15V	FCP-9P	
	FG-24V	FB-18V	FI-7BS	
	FG-29V	FB-19TV	FI-9BS	

### OF-D6 Eyepiece Cap

<b>Model Name</b> OF-D6	<b>Description</b> Eyepiece cap		
Compatible Endoscopes	FB-15RBS	FI-9RBS	FI-13RBS
	FB-18RBS	FI-10RBS	FI-16RBS

### OF-D1 Eyepiece Cap



Model Name	Description			
OF-D1	Eyepiece cap			
Compatible Endoscopes	FD-34V2	FB-8V	FB-19TV	FI-13BS
	FG-16V	FB-10V	FCN-15X	FI-16BS
	FG-24V	FB-15V	FI-9BS	
	FG-29V	FB-18V	FI-10BS	

### OF-D3 Eyeshield



Model Name	Description			
OF-D3	Eyeshield			
Compatible Endoscopes	FD-34V2	FG-24V	FB-15V	FB-19TV
	FG-16V	FG-29V	FB-18V	FCN-15X

### OF-C3 Soaking Cap



Model Name	Description
OF-C3	Soaking cap
Compatible Endoscopes	FCN-15X

### OF-B82 Adjustable Luer Lock Seal



Model Name	Description	
OF-B82	Adjustable luer lock seal	
Compatible Endoscopes	FHY-10RBS	FHY-15RBS

**GM00261 Air/Water Valve**

<b>Model Name</b> GM00261	<b>Description</b> Air/water valve for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

**GM00259 Suction Valve**

<b>Model Name</b> GM00259	<b>Description</b> Suction valve for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

**MF00004 Biopsy Valve**

<b>Model Name</b> MF00004	<b>Description</b> Biopsy valve for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

**GS00006 Water Bottle**

<b>Model Name</b> GS00006	<b>Description</b> Water bottle for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL



#### Manufacturers

##### VERSA Consumables

 **PENTAX-AOHUA Medical Technologies Co., Ltd.**  
East of 3rd Floor, Block C, Building 1,  
No.5 Shenwang Road,  
Minhang District, Shanghai 201108,  
P. R. China.

##### All other products

 **HOYA Corporation**  
6-10-1, Nishi-shinjuku  
Shinjuku-ku, Tokyo  
160-0023 Japan

 **EMEA Headquarter  
Germany**

PENTAX Europe GmbH  
Julius Vosseler Straße 104, 22527 Hamburg, Germany



Medical Devices Class: I

 **EMEA Headquarter  
Germany**

PENTAX Europe GmbH  
Julius Vosseler Straße 104, 22527 Hamburg, Germany



0123  
TÜV SÜD CE 0123 Medical Devices Class: IIa, IIb, Is and I

#### Distributor

##### EMEA Headquarter Germany

PENTAX Europe GmbH  
Julius Vosseler Straße 104  
22527 Hamburg, Germany  
Tel.: +49 40 / 5 61 92 - 0  
Fax: +49 40 / 5 60 42 13  
E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative

**PENTAX**  
**MEDICAL**