

Hygiene Consumables

06

Hygiene Consumables

table of contents

06-1 PROfile Single-Use Cleaning Brushes

06-2 Cleaning Brushes

06-3 Reprocessing Accessories

06-4 VERSA Consumables



DESIGNED FOR PENTAX Medical – One product range that fits every PENTAX Medical endoscope

Manually pre-cleaning flexible endoscopes is key to ensuring reliable reprocessing as an end result. Even with magnification, visual assessment of the manual pre-cleaning results is insufficient to detect microorganisms. As such, PENTAX Medical provides a verified and validated solution to improve confidence in manual pre-cleaning, with the new PROfILE single-use cleaning brushes.

1

Proven cleaning efficiency

2








Designed and validated to meet PENTAX Medical's quality standards

3



Compliant with the industry's highest reprocessing standards

Overview PROfile single-use cleaning brushes

Long instrument/suction channel brushes



Short instrument channel inlet brushes




Find the right brush with the brush finder on www.the-hygiene-solution-that-fits.com

Proven cleaning efficiency

In combination with PENTAX Medical endoscopes, the PROfile single-use cleaning brushes were tested against residual soil and compared against the performance of legacy brushes.

Through an authorized statistical test method, the internal validation proved both brushes to have a similar performance.



Designed and validated to meet PENTAX Medical's quality standards

PROfile single-use cleaning brushes were made to meet the high-quality standards of PENTAX Medical engineering. Extensive internal tests have proven that they are the best possible fit to ensure proper hygiene results, with no impairment to the endoscope, across all types.

- Tailored for PENTAX Medical endoscopes' precise dimensions, structures, and materials
- Improved bonding and pull-off capability at optimal insertion resistance of $\leq 10N$
- No damage to endoscope working channels and reliable even at increased thermal conditions

Compliant with the industry's highest reprocessing standards

How clean is clean enough? This was the leading question of PENTAX Medical engineers when deciding upon the validation criteria for hygienic efficacy of the PROfile single-use cleaning brushes. In order to follow established standards, PENTAX Medical validated against the requirements of the Association for the Advancement of Medical Instrumentation (AAMI)¹.




Article	Hemoglobin [< 2.2µg/cm ²]			Protein [< 6.4µg/cm ²]			Carbohydrate [< 1.8µg/cm ²]		
	µg/article	µg/cm ²	Result	µg/article	µg/cm ²	Result	µg/article	µg/cm ²	Result
Endoscope 1	14.88	1.06	PASS	31.05	2.20	PASS	10.82	0.77	PASS
Endoscope 2	9.92	0.70	PASS	22.68	1.61	PASS	8.05	0.57	PASS
Endoscope 3	19.84	1.41	PASS	59.40	4.21	PASS	6.82	0.48	PASS

Test example: Residual of three markers at scope's distal end after third cleaning cycle, complying with AAMI requirements.

PROfILE Single-Use Cleaning Brushes

One product range that matches every channel

Brush				Compatible endoscope							
Model	Length [mm]	Outer Diameter [ø mm]	Color [Shaft / Handle]	Compatible Channel I.D. [mm]	Working Length [mm]	Gastro-ent-erology	Pulmo-nology	Urology	ENT & Speech	Intu-bation	EUS
CS2010A	1,020	2	Blue	1.1 - 1.5	Up to 720		•	•	•	•	
CS2022A	2,200	2		1.1 - 1.5	Up to 1,900						•
CS2037A	3,720	2		1.1 - 1.5	Up to 3,420	•					
CS3022A	2,220	3	White	1.7 - 2.4	Up to 1,920	•					
CS4010A	1,020	4	Orange	1.8 - 3.2	Up to 720		•	•	•	•	
CS5503A	310	5.5	Yellow	2.0 - 5.0	N/A		•	•	•	•	
CS5522A	2,420	5.5		2.0 - 5.0	Up to 2,120	•					
CS5530A	3,020	5.5		2.0 - 5.0	Up to 2,720	•					
CS-C11A	N/A	6-14	Cobalt Blue	N/A	N/A		•	•	•	•	
CS-C12A	N/A	5	Yellow	N/A	N/A						•
CS-C13A	N/A	6-15	Green	N/A	N/A	•					



PROFILE Single-Use Cleaning Brushes

One product range to match every channel



PENTAX Medical PROFILE single-use endoscope cleaning brush kits with proven usability are available in a wide range of dimensions to cover all channel sizes.

Compliant with reprocessing guidelines

- Single-use cleaning accessories are internationally recommended by reprocessing guidelines.

Meeting infection control requirements

- PENTAX Medical single-use cleaning brush kits help to ensure safe reprocessing and simplify processes.

Proven cleaning results

- Streamlined product range to match all PENTAX Medical endoscopes in all product segments
- Tailored for PENTAX Medical endoscopes precise dimensions, structures and materials
- Validated cleaning efficacy

PROFILE Shaft Color	Model Name	Working Length [mm]	Endoscope Channel Diameter [mm]	Endoscope type (both video and fiberoptic)					
				Gastroen- terology	Pulmono- logy	Urology	ENT & Speech	Intuba- tion	EUS
Blue	CS2010A	Up to 720	1.1 to 1.5		•	•	•	•	
	CS2022A	Up to 1,900	1.1 to 1.5						•
	CS2037A	Up to 3,420	1.1 to 1.5	•					
White	CS3022A	Up to 1,920	1.7 to 2.4	•					
Orange	CS4010A	Up to 720	1.8 to 3.2		•	•	•	•	
Yellow	CS5503A	N/A	2.0 to 5.0		•	•	•	•	
	CS5522A	Up to 2,120	2.0 to 5.0	•					
	CS5530A	Up to 2,720	2.0 to 5.0	•					
Cobalt Blue	CS-C11A	N/A	N/A		•	•	•	•	
Yellow	CS-C12A	N/A	N/A						•
Green	CS-C13A	N/A	N/A	•					

Consult your local representative for information on the products currently available in your region.

PENTAX Medical

Accessories Programme for cleaning

for your PENTAX Medical endoscopes

Original is better

Finely tuned system components are essential for safe and successful endoscopy. That is why PENTAX Medical original accessories aim to provide the highest possible quality in every detail. The following pages describe our carefully compiled accessory programme to ensure safe and effective cleaning of your PENTAX Medical endoscope.



How to choose the correct accessories:

Before you decide on new accessories, please ensure that your PENTAX Medical endoscope is compatible with the desired accessory.

- The illustration below is intended to help you to identify the relevant areas of your endoscope quickly and simply. Our PENTAX Medical specialists will be pleased to advise you if you have any additional questions.
- All PENTAX Medical original accessories in this catalogue have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

Orientation

In order to make your choice of PENTAX Medical accessories as simple and reliable as possible, we have identified fundamental product features with simple-to-understand symbols displayed next to the relevant technical data.




Autoclavable



**For single
use only**




CS6021T Tri-Bristled Cleaning Brushes (10 pcs)



Model name	Description			
CS6021T	Single-use cleaning brushes for channel, 10 pcs.			
Compatible endoscopes	EC-2990Li	EC-3890Mi2	EC34-i10TM	EG-2990Zi
	EC-2990Fi	EC-3890LK	EC38-i10L	EG-3490K
	EC-2990Mi	EC-3890FK	EC38-i10F	EG-3890TK
	EC-3490Li	EC-3890MK	EC38-i10M	EG16-K10
	EC-3490Fi	EC-3890FK2	EC38-i10F2	EG27-i10
	EC-3490Mi	EC-3890MK2	EC38-i10M2	EG29-i10
	EC-3490LK	EC-3890LZi	EC38-i10NL	EG34-i10
	EC-3490FK	EC-3890FZi	EC38-i10NF	EE-1580K
	EC-3490MK	EC-3890MZi	ED-3490TK	EG-1870K
	EC-3490TLi	EC-3890TLK	ED34-i10T	EG-3270UK
	EC-3490TFi	EC-3890TFK	EG-1690K	EG-3670URK
	EC-3490TMi	EC34-i10L	EG-2490K	EG-3870UTK
	EC-3890Li	EC34-i10F	EG-2790i	ES-3870K
	EC-3890Fi	EC34-i10M	EG-2790K	
	EC-3890Mi	EC34-i10TL	EG-2990i	
	EC-3890Fi2	EC34-i10TF	EG-2990K	




CS3010S Cleaning Brush for Endoscope Channel



Model name	Description			
CS3010S	Cleaning brush for channel			
Compatible endoscopes	EB-1170K	FB-10V	FI-9RBS	FI-13BS
	EB-1990i	FHY-10RBS	FI-10BS	FI-13RBS
	FB-8V	FI-9BS	FI-10RBS	FUR-9RBS



CS3025S Cleaning Brush for Endoscope Channel



Model name	Description	
CS3025S	Cleaning brush for A/W and WJ channel	
Compatible endoscopes	EE-1580K	FCP-9P




CS6002SN Cleaning Brush for Endoscope Channel



Model name	Description			
CS6002SN	Cleaning brush for channel			
Compatible endoscopes	EB-1170K	EB-1990i	FB-15RBS	FI-10BS
	EB-1570AK	EB-1970UK	FB-18BS	FI-10RBS
	EB-1570K	ECY-1570K	FB-18RBS	FI-13BS
	EB-1572K	ECY-1575K	FCN-15X	FI-13RBS
	EB-1575K	FB-8V	FCP-9P	FI-16BS
	EB-1970AK	FB-10V	FCY-15RBS	FI-16RBS
	EB-1970K	FB-15V	FHY-10RBS	FNL-15RP3
	EB-1970TK	FB-18V	FHY-15RBS	FUR-9RBS
	EB-1972K	FB-19TV	FI-9BS	VNL-1570STK
	EB-1975K	FB-15BS	FI-9RBS	




CS6015ST Cleaning Brush for Endoscope Channel



Model name	Description			
CS6015ST	Cleaning brush for channel			
Compatible endoscopes	FG-16V	EB-1970K	FB-15V	FCN-15X
	FG-24V	EB-1970TK	FB-18V	FCY-15RBS
	EB-1570AK	EB-1972K	FB-19TV	FHY-15RBS
	EB-1570K	EB-1975K	FB-15BS	FI-16BS
	EB-1572K	EB-1970UK	FB-15RBS	FI-16RBS
	EB-1575K	ECY-1570K	FB-18BS	FNL-15RP3
	EB-1970AK	ECY-1575K	FB-18RBS	VNL-1570STK




CS6021SN Cleaning Brush for Endoscope Channel



Model name	Description			
CS6021SN	Cleaning brush for working channel			
Compatible endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	FG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V



CS6030SN Cleaning Brush for Endoscope Channel



Model name CS6030SN	Description Cleaning brush for working channel, extra long
Compatible endoscopes	VSB-2990i



CS-C3S Cleaning Brush for Endoscope Cylinder



Model name CS-C3S	Description Cleaning brush for cylinders or certain valves			
Compatible endoscopes	EB-1170K	EB-1990i	FB-15RBS	FI-10BS
	EB-1570AK	EB-1970UK	FB-18BS	FI-10RBS
	EB-1570K	ECY-1570K	FB-18RBS	FI-13BS
	EB-1572K	ECY-1575K	FCN-15X	FI-13RBS
	EB-1575K	FB-8V	FCP-9P	FI-16BS
	EB-1970AK	FB-10V	FCY-15RBS	FI-16RBS
	EB-1970K	FB-15V	FHY-10RBS	FNL-15RP3
	EB-1970TK	FB-18V	FHY-15RBS	FUR-9RBS
	EB-1972K	FB-19TV	FI-9BS	VNL-1570STK
	EB-1975K	FB-15BS	FI-9RBS	



CS-C5S Cleaning Brush for Endoscope Cylinder



Model name CS-C5S	Description Cleaning brush for cylinder			
Compatible endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	FG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V



CS-C9S Cleaning Brush for Endoscope Cylinder



Model name CS-C9S	Description Cleaning brush for valve cylinders, inlet/channel ports, scope tips/ elevator mechanisms and valves			
Compatible endoscopes	EC-2990Li	EC-3890Mi2	EC34-i10TM	EG-2990Zi
	EC-2990Fi	EC-3890LK	EC38-i10L	EG-3490K
	EC-2990Mi	EC-3890FK	EC38-i10F	EG-3890TK
	EC-3490Li	EC-3890MK	EC38-i10M	EG16-K10
	EC-3490Fi	EC-3890FK2	EC38-i10F2	EG27-i10
	EC-3490Mi	EC-3890MK2	EC38-i10M2	EG29-i10
	EC-3490LK	EC-3890LZi	EC38-i10NL	EG34-i10
	EC-3490FK	EC-3890FZi	EC38-i10NF	VSB-2990i
	EC-3490MK	EC-3890MZi	ED-3490TK	EE-1580K
	EC-3490TLi	EC-3890TLK	ED34-i10T	EG-1870K
	EC-3490TFi	EC-3890TFK	EG-1690K	EG-3270UK
	EC-3490TMi	EC34-i10L	EG-2490K	EG-3670URK
	EC-3890Li	EC34-i10F	EG-2790i	EG-3870UTK
	EC-3890Fi	EC34-i10M	EG-2790K	ES-3870K
	EC-3890Mi	EC34-i10TL	EG-2990i	
	EC-3890Fi2	EC34-i10TF	EG-2990K	



OF-B212 Cleaning Adapter



Model name	Description		
OF-B212	Cleaning adapter		
Compatible endoscopes	EB11-J10	EB19-J10	VNL15-J10
	EB15-J10	EB19-J10U	

Fitting O-ring set is OF-B214



OF-G17 Air-Water Channel Cleaning Adapter



Model name	Description			
OF-G17	A/W cleaning adapter			
Compatible endoscopes	EC-2990Li	EC-3890FK	EC38-i10M	EG-3490K
	EC-2990Fi	EC-3890MK	EC38-i10F2	EG-3890TK
	EC-2990Mi	EC-3890FK2	EC38-i10M2	EG16-K10
	EC-3490Li	EC-3890MK2	EC38-i10NL	EG17-J10
	EC-3490Fi	EC-3890LZi	EC38-i10NF	EG27-i10
	EC-3490Mi	EC-3890FZi	EC38-i10cL	EG29-i10
	EC-3490LK	EC-3890MZi	EC38-i10cF	EG29-i10c
	EC-3490FK	EC-3890TLK	EC38-i10cM	EG34-i10
	EC-3490MK	EC-3890TFK	EC38-i10cF2	FC-38LV
	EC-3490TLi	EC34-i10L	ED-3490TK	FC-38FV
	EC-3490TFi	EC34-i10F	ED32-i10	FD-34V2
	EC-3490TMi	EC34-i10M	ED34-i10T	FG-16V
	EC-380LKp	EC34-i10NF	ED34-i10T2	FG-24V
	EC-380FKp	EC34-i10NL	EE17-J10	FG-29V
	EC-380FK2p	EC34-i10TL	EG-1690K	FS-34V
	EC-380MK2p	EC34-i10TF	EG-2490K	VSB-2990i
	EC-3890Li	EC34-i10TM	EG-2790i	EG-1870K
	EC-3890Fi	EC34-i10cL	EG-2790K	EG-3270UK
	EC-3890Mi	EC34-i10cF	EG-290Kp	EG-3870UTK
	EC-3890Fi2	EC34-i10cM	EG-2990i	ES-3870K
	EC-3890Mi2	EC38-i10L	EG-2990K	
	EC-3890LK	EC38-i10F	EG-2990Zi	



OE-C20 Water Jet Channel Cleaning Adapter



Model name	Description			
OE-C20	Water jet channel cleaning adapter			
Compatible endoscopes	EC-2990Li	EC-3890Fi2	EC34-i10NF	EC38-i10cL
	EC-2990Fi	EC-3890Mi2	EC34-i10NL	EC38-i10cF
	EC-2990Mi	EC-3890LK	EC34-i10TL	EC38-i10cM
	EC-3490Li	EC-3890FK	EC34-i10TF	EC38-i10cF2
	EC-3490Fi	EC-3890MK	EC34-i10TM	EG-2990i
	EC-3490Mi	EC-3890FK2	EC34-i10cL	EG-2990K
	EC-3490LK	EC-3890MK2	EC34-i10cF	EG-2990Zi
	EC-3490FK	EC-3890LZi	EC34-i10cM	EG-3490K
	EC-3490MK	EC-3890FZi	EC38-i10L	EG-3890TK
	EC-3490TLi	EC-3890MZi	EC38-i10F	EG29-i10
	EC-3490TFi	EC-3890TLK	EC38-i10M	EG29-i10c
	EC-3490TMi	EC-3890TFK	EC38-i10F2	EG34-i10
	EC-3890Li	EC34-i10L	EC38-i10M2	VSB-2990i
	EC-3890Fi	EC34-i10F	EC38-i10NL	ES-3870K
	EC-3890Mi	EC34-i10M	EC38-i10NF	



OF-B153 Air-Water Suction Channel Cleaning Adapter



Model name	Description			
OF-B153	Cleaning adapter			
Compatible endoscopes	EC-2990Li	EC-3890FK	EC34-i10cF	EG-1690K
	EC-2990Fi	EC-3890MK	EC34-i10cM	EG-2490K
	EC-2990Mi	EC-3890FK2	EC38-i10L	EG-2790i
	EC-3490Li	EC-3890MK2	EC38-i10F	EG-2790K
	EC-3490Fi	EC-3890LZi	EC38-i10M	EG-2990i
	EC-3490Mi	EC-3890FZi	EC38-i10F2	EG-2990K
	EC-3490LK	EC-3890MZi	EC38-i10M2	EG-2990Zi
	EC-3490FK	EC-3890TLK	EC38-i10NL	EG-3490K
	EC-3490MK	EC-3890TFK	EC38-i10NF	EG-3890TK
	EC-3490TLi	EC34-i10L	EC38-i10cL	EG16-K10
	EC-3490TFi	EC34-i10F	EC38-i10cF	EG17-J10
	EC-3490TMi	EC34-i10M	EC38-i10cM	EG27-i10
	EC-3890Li	EC34-i10NF	EC38-i10cF2	EG29-i10
	EC-3890Fi	EC34-i10NL	ED-3490TK	EG29-i10c
	EC-3890Mi	EC34-i10TL	ED32-i10	EG34-i10
	EC-3890Fi2	EC34-i10TF	ED34-i10T	VSB-2990i
	EC-3890Mi2	EC34-i10TM	ED34-i10T2	EG-3670URK
	EC-3890LK	EC34-i10cL	EE17-J10	



OF-B218 Cleaning Cap for Instrument Channel Inlet



Model name OF-B218	Description Cleaning cap for instrument channel inlet			
Compatible endoscopes	EB11-J10	EB19-J10	EG34-J10U	EG38-J10UT
	EB15-J10	EB19-J10U	EG36-J10UR	VNL15-J10



OE-B11 PENTAX Medical Cleaning Adapter



Model name OE-B11	Description Cleaning adapter			
Compatible endoscopes	EB19-J10U			



OE-B10 Cleaning Adapter



Model name OE-B10	Description Cleaning adapter			
Compatible endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT	

Fitting O-ring set is OF-B195



OE-G1 Cleaning Adapter



Model name OE-G1	Description Cleaning adapter			
Compatible endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT	

OE-C9 PVE Soaking Cap



Model name	Description			
OE-C9	Soaking cap			
Compatible endoscopes	EC-2990Li	EC-3890FZi	EE17-J10	EB-1990i
	EC-2990Fi	EC-3890MZi	EG-1690K	EB11-J10
	EC-2990Mi	EC-3890TLK	EG-2490K	EB15-J10
	EC-3490Li	EC-3890TFK	EG-2790i	EB19-J10
	EC-3490Fi	EC34-i10L	EG-2790K	EB-1970UK
	EC-3490Mi	EC34-i10F	EG-290Kp	EB19-J10U
	EC-3490LK	EC34-i10M	EG-2990i	ECY-1575K
	EC-3490FK	EC34-i10NF	EG-2990K	EE-1580K
	EC-3490MK	EC34-i10NL	EG-2990Zi	EG-1870K
	EC-3490TLi	EC34-i10TL	EG-3490K	EG-3270UK
	EC-3490TFi	EC34-i10TF	EG-3890TK	EG-3670URK
	EC-3490TMi	EC34-i10TM	EG16-K10	EG-3870UTK
	EC-380LKp	EC38-i10L	EG17-J10	EG34-J10U
	EC-380FKp	EC38-i10F	EG27-i10	EG36-J10UR
	EC-380FK2p	EC38-i10M	EG29-i10	EG38-J10UT
	EC-380MK2p	EC38-i10F2	EG34-i10	ES-3870K
	EC-3890Li	EC38-i10M2	VSF-2990i	VNL-1070STK
	EC-3890Fi	EC38-i10NL	EB-1170K	VNL-1190STK
	EC-3890Mi	EC38-i10NF	EB-1570AK	VNL-1570STK
	EC-3890Fi2	EC38-i10cL	EB-1570K	VNL-1590STi
	EC-3890Mi2	EC38-i10cF	EB-1572K	VNL8-J10
	EC-3890LK	EC38-i10cM	EB-1575K	VNL11-J10
	EC-3890FK	EC38-i10cF2	EB-1970AK	VNL15-J10
	EC-3890MK	ED-3490TK	EB-1970K	PVK-1070Z
	EC-3890FK2	ED32-i10	EB-1970TK	PVK-J10
	EC-3890MK2	ED34-i10T	EB-1972K	
	EC-3890LZi	ED34-i10T2	EB-1975K	

OE-C27 Soaking Cap for Navigation Connector



Model name	Description			
OE-C27	Soaking cap for navigation connector			
Compatible endoscopes	EC34-i10NF	EC34-i10NL	EC38-i10NF	EC38-i10NL

OE-U1 Soaking Cap for Scanning Unit Connector



Model name	Description			
OE-U1	Soaking cap for scanning unit connector			
Compatible endoscopes	EB-1970UK	EG-3270UK	EG-3670URK	EG-3870UTK

OE-U4 Soaking Cap for PENTAX Medical Ultrasound Endoscopes



Model name	Description			
OE-U4	Soaking cap for PENTAX Medical ultrasound endoscopes			
Compatible endoscopes	EB19-J10U	EG34-J10U	EG36-J10UR	EG38-J10UT



OF-B162 Cleaning Cap for 2ch



Model name	Description	
OF-B162	Cleaning cap for 2ch	
Compatible endoscopes	EC-3890TLK	EC-3890TFK



OF-B115 Suction & Air-Water Channel Cleaning Adapter



Model name	Description			
OF-B115	Cleaning cap for A/W/suction channel			
Compatible endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	EG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V



OF-B166 Air-Water-Suction Channel Cleaning Adapter



Model name OF-B166	Description Cleaning adapter		
Compatible endoscopes	EE-1580K	EG-1870K	ES-3870K



OE-C8 Suction Nipple Cleaning Adapter



Model name OE-C8	Description Cleaning adapter suction nipple		
Compatible endoscopes	EG-1870K	ES-3870K	



OF-B170 Air Cleaning Brush Adapter



Model name OF-B170	Description Cleaning brush adapter for A/W channel EG-1580K		
Compatible endoscopes	EE-1580K		



OE-C5 Air-Water & Suction Channel Adapter



Model name OE-C5	Description Cleaning adapter for A/W channel		
Compatible endoscopes	EE-1580K	EG-3670URK	



OF-B103 Suction Channel Cleaning & Disinfection Adapter



Model name	Description			
OF-B103	Suction channel cleaning & disinfection adapter			
Compatible endoscopes	FB-15BS	FCY-15RBS	FI-10BS	FI-16RBS
	FB-15RBS	FHY-10RBS	FI-10RBS	FNL-15RP3
	FB-18BS	FHY-15RBS	FI-13BS	FUR-9RBS
	FB-18RBS	FI-9BS	FI-13RBS	
	FCP-9P	FI-9RBS	FI-16BS	



BS-B1 Soaking Cap for BS



Model name	Description			
BS-B1	Soaking cap for BS series			
Compatible endoscopes	FB-15BS	FHY-10RBS	FI-9RBS	FI-16BS
	FB-15RBS	FHY-15RBS	FI-10BS	FI-16RBS
	FB-18BS	FI-7BS	FI-10RBS	FNL-10RBS
	FB-18RBS	FI-7RBS	FI-13BS	FUR-9RBS
	FCY-15RBS	FI-9BS	FI-13RBS	



OF-B155 Cleaning Adapter



Model name	Description			
OF-B155	Cleaning cap for suction channel			
Compatible endoscopes	EB-1170K	EB-1575K	EB-1972K	ECY-1575K
	EB-1570AK	EB-1970AK	EB-1975K	VNL-1570STK
	EB-1570K	EB-1970K	EB-1990i	
	EB-1572K	EB-1970TK	EB-1970UK	



OF-B75 Suction Channel Cleaning Adapter



Model name	Description		
OF-B75	Cleaning cap for suction channel		
Compatible endoscopes	FB-8V	FB-15V	FB-19TV
	FB-10V	FB-18V	FCN-15X

SHA-P6 Leakage Tester



Model name
SHA-P6

Description
Leakage tester

GS00003 Cleaning Kit

Model name	Description			
GS00003	Cleaning kit for VERSA endoscopes			
Compatible endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

GF00089 Cleaning Plug

Model name	Description			
GF00089	Cleaning plug for VERSA endoscopes			
Compatible endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

GF00027 Leakage Tester

Model name	Description			
GF00027	Leakage tester for VERSA endoscopes			
Compatible endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

GM00906 Water Resistant Cap


Model Name	Description			
GM00906	Water resistant cap for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL


GS00017 Plug Label

Model Name	Description			
GS00017	Plug label for VERSA endoscopes			
Compatible Endoscopes	EG27-V10c	EG29-V10c	EC38-V10cM	EC38-V10cL

Manufacturers

VERSA Consumables

 **PENTAX-AOHUA Medical Technologies Co., Ltd.**
East of 3rd Floor, Block C, Building 1,
No.5 Shenwang Road,
Minhang District, Shanghai 201108,
P. R. China.

 **EMEA Headquarter
Germany**

PENTAX Europe GmbH
Julius Vosseler Straße 104, 22527 Hamburg, Germany



Medical Devices Class: I

All other products

 **HOYA Corporation**
6-10-1, Nishi-shinjuku
Shinjuku-ku, Tokyo
160-0023 Japan

 **EMEA Headquarter
Germany**

PENTAX Europe GmbH
Julius Vosseler Straße 104, 22527 Hamburg, Germany



0123
TÜV SÜD CE 0123 Medical Devices Class: I, Is, IIa and IIb

Distributor

EMEA Headquarter Germany

PENTAX Europe GmbH
Julius Vosseler Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL