


Therapeutic Accessories – The Choice for Quality and Efficiency in Endoscopy.

Available to ensure highest standards.



PENTAX Medical

State-of-the-art accessories product line for diagnosis and treatment in endoscopy.



Our commitment to quality

Safe and successful endoscopy depends upon the use of a finely tuned endoscopic system of components. This is the underlying principle behind PENTAX Medical's commitment to a product line of therapeutic accessories with the highest possible level of quality. We invite you to examine our carefully selected and meticulously crafted therapeutic accessories for state-of-the-art diagnostic and treatment applications. All PENTAX Medical accessories have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

How to select the correct accessories.

Before you select new accessories, please consider the following:

- The compatibility of your endoscope with the desired diagnostic/treatment application.
- The diameter of the instrument channel and working length of your endoscope.
- Use only accessories that are compatible with your endoscope to ensure the highest standards of endoscopic therapy.

PENTAX Medical

Single Use Accessories


Single Use

Biopsy Forceps




Description	Model Name	Compatible Endoscopes	
		Instrument Channel Diameter [mm]	Working Length [mm]
with Windows	KW-D1810	2.0	1000
	KW-D1816	2.0	1600
	KW-D2416T	2.8	1600
	KW-D2423T	2.8	2300
with Windows/ Spike	KH-D2416T	2.8	1600
	KH-D2423T	2.8	2300
Crocodile	KA-D2416T	2.8	1600
	KA-D2423T	2.8	2300
Crocodile with Spike	KB-D2416T	2.8	1600
	KB-D2423T	2.8	2300

Single Use Retrieval Basket



Model Name	Compatible Endoscopes		Width of Basket [mm]
	Instrument Channel Diameter [mm]	Working Length [mm]	
GB-D1807L015	2.0	700	15
GB-D1810L020	2.0	1000	20
GB-D1819L020	2.0	1900	20
GB-D2415L030	2.8	1500	30
GB-D2419L020	2.8	1900	20
GB-D2423L035	2.8	2300	35

Single Use Polypectomy Snare



Description	Model Name	Compatible Endoscopes		Width of Snare Wire [mm]
		Instrument Channel Diameter [mm]	Working Length [mm]	
Oval	DO-D2416-15	2.8	1650	15
	DO-D2423-15	2.8	2300	15
	DO-D2416-20	2.8	1650	20
	DO-D2423-20	2.8	2300	20
	DO-D2416-25	2.8	1650	25
	DO-D2423-25	2.8	2300	25
Hexagonal	DH-D2416-15	2.8	1650	15
	DH-D2423-15	2.8	2300	15
	DH-D2423-20	2.8	2300	20
	DH-D2416-25	2.8	1650	25
	DH-D2423-25	2.8	2300	25


Single Use Injection Needle



Needle Size	Model Name	Compatible Endoscopes		Needle Length [mm]
		Instrument Channel Diameter [mm]	Working Length [mm]	
23G	NI-D1816-T2304	2.0	1650	4
	NI-D1816-T2305	2.0	1650	5
	NI-D1816-T2306	2.0	1650	6
	NI-D1816-T2308	2.0	1650	8
	NI-D2416-T2304	2.8	1650	4
	NI-D2416-T2306	2.8	1650	6
	NI-D2423-T2304	2.8	2300	4
	NI-D2423-T2305	2.8	2300	5
	NI-D2423-T2306	2.8	2300	6
25G	NI-D1816-T2504	2.0	1650	4
	NI-D1816-T2505	2.0	1650	5
	NI-D2416-T2504	2.8	1650	4
	NI-D2416-T2505	2.8	1650	5
	NI-D2416-T2506	2.8	1650	6
	NI-D2423-T2504	2.8	2300	4
	NI-D2423-T2505	2.8	2300	5
	NI-D2423-T2506	2.8	2300	6

Single Use

Sprayer



Model Name	Compatible Endoscopes	
	Instrument Channel Diameter [mm]	Working Length [mm]
TJ-D2418PB	2.8	1800
TJ-D2423PB	2.8	2300

PENTAX Medical

Reusable Accessories (Autoclavable)

Reusable

Crocodile Biopsy Forceps



Model Name	Compatible Endoscopes		Note
	Instrument Channel Diameter [mm]	Working Length [mm]	
KA1806S	1.8	600	Autoclavable
KA1811S	2.0	1100	
KA1815S	2.0	1500	
KA2411S	2.8	1100	
KA2415S	2.8	1500	
KA2422S	2.8	2200	

Reusable

Biopsy Forceps with Window



Model Name	Compatible Endoscopes		Note	
	Instrument Channel Diameter [mm]	Working Length [mm]		
KW1511S	2.0	1100	Autoclavable	
KW1806S	2.0	600		
KW1811S	2.0	1100		
KW1815S	2.0	1500		
KW1818S	2.0	1800		
KW1818CS	2.0	1800		for ERCP scopes
KW2211S	2.4	1100		
KW2215S	2.4	1500		
KW2218CS	2.4	1800		for ERCP scopes
KW2411S	2.8	1100		
KW2415S	2.8	1500		
KW2418S	2.8	1800		
KW2422S	2.8	2200		
KW2430S	2.8	3000		

Reusable**Biopsy Forceps with Window and Spike**

Model Name	Compatible Endoscopes		Note	
	Instrument Channel Diameter [mm]	Working Length [mm]		
KH2411S	2.8	1100	Autoclavable	
KH2415S	2.8	1500		
KH2418CS	2.8	1800		for ERCP scopes
KH2422S	2.8	2200		
KH2430S	2.8	3000		

Reusable**Biopsy Forceps with Mini Cups**

Model Name	Compatible Endoscopes		Note
	Instrument Channel Diameter [mm]	Working Length [mm]	
KS1011CS	1.2	1100	Autoclavable
KS1022CS	1.2	2200	

Reusable**Biopsy Forceps with Window and Extra Large Cups**

Model Name	Compatible Endoscopes		Note
	Instrument Channel Diameter [mm]	Working Length [mm]	
KW3111S	3.2	1100	Autoclavable
KW3115S	3.2	1500	
KW3122S	3.2	2200	

Reusable**Biopsy Forceps with Round Cups**

Model Name	Compatible Endoscopes		Note
	Instrument Channel Diameter [mm]	Working Length [mm]	
KW2411R	2.8	1100	Autoclavable
KW2415R	2.8	1500	
KW2418R	2.8	1800	
KW2422R	2.8	2200	
KW2430R	2.8	3000	

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg
Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: medical@pentax.de

PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
United Kingdom
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical@pentax.co.uk

PENTAX France Life Care S.A.S.
112 quai de Bezons
B.P. 204
95106 ARGENTEUIL CEDEX
France
Tel.: +33 1 / 30 25 75 75
Fax: +33 1 / 30 25 74 45
E-mail: contact.medical@pentax.fr

PENTAX Nederland B.V.
Amsterdamseweg 29
1420 CA Uithoorn
Netherlands
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: medical@pentax.nl

PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Italy
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.lifecare@pentaxitalia.it

**SIMMEDICA – Sistemas Integrales
de Medicina, S.A.**
Avenida del Sistema Solar 25
28830 San Fernando de Henares · Madrid
Spain
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

PENTAX Europe GmbH
Representative office in Moscow
Sadovnicheskaya str, 82, build 2, entrance 6
Regus, Business center „Aurora“
Office 2012, 2013
Moscow, 115035
Russian Federation
Tel.: +7 495 / 792 52 00
Fax: +7 495 / 792 35 66

**PENTAX Europe GmbH
Turkey Liaison Office**
Veko Giz Plaza, Meydan Sok. No:3/43
34396 Maslak – Istanbul
Turkey
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

**HOYA Corporation
PENTAX Life Care Division**
1-1-110, Tsutsujigaoka
Akishima-shi
196-0012 Tokyo
Japan
Tel.: +81 3 / 39 60 51 55
Fax: +81 3 / 53 92 67 24